
Mediflo Duo

EN - RESPIRATORY EXERCISER FOR INSPIRATION EXERCISES (SMI) AND EXPIRATION EXERCISES (PEP)

Intended Purpose

Breathing Trainers are manually operated incentive spirometers intended for inhalation training according to the Sustained Maximal Inspiration (SMI) method and expiration training according to the Positive Expiratory Pressure (PEP) method.

They are intended for use in adult and paediatric patients by medical professionals, carers and patients in hospital and out-of-hospital settings. Breathing Trainers are nonsterile single patient use devices, to be used for up to 30 days in a hospital environment and up to 6 months in a home environment.

Indications

Indications for the use of Breathing Trainers include:

- Acute and chronic respiratory failure caused by surgery
- Neurological or musculoskeletal dysfunction, e.g., MND
- Increasing lung volume by increasing FRC and Vt
- Reducing hyperinflation / air trapping in e.g.: emphysema, bronchitis, asthma
- Improving airway clearance in; cystic fibrosis (over 4-year-old); chronic bronchitis; bronchiectasis; bronchiolitis obliterans
- Maximising the delivery of bronchodilators in patients receiving bronchial hygiene therapy

Contraindications

Contraindications for the use of Breathing Trainers include:

- Untreated pneumothorax
- Intracranial pressure > 20mm Hg
- Active haemoptysis
- Recent trauma or surgery to the skull, face, mouth, or oesophagus
- Patient with acute asthma attack or acute worsening of Chronic Obstructive Pulmonary Disease (COPD), unable to tolerate increased work of breathing
- Acute sinusitis or epistaxis
- Tympanic membrane rupture or other known or suspected inner ear pathology
- Nausea
- Following coronary artery bypass graft or valve replacement surgery
- When intended user:
 - is unable to generate adequate inspiration with a vital capacity < 10 mL/kg or an inspiratory capacity < 33% of predicted normal range.
 - unable to deep breathe effectively due to pain, diaphragmatic dysfunction, or opiate analgesia
 - cannot be instructed or supervised to assure appropriate use of the device, or where co-operation and user understanding of the device is absent.
 - is confused, delirious, heavily sedated or comatose.

Precautions/Cautions

- Visually inspect the product prior to use and ensure that all connections are secure.
- Do not use if there is damage, discoloration or anything that may reduce the functionality of the device.
- Do not use if package has been previously opened or damaged.
- Do not use after expiry date.

Warnings

- Frequency of use and recommended inspiratory/ expiratory resistance should be performed at the direction of doctor.
- Consult doctor before use if the patient is on medication or suffer from a medical condition. Do not stop taking any medication unless advised by doctor.

- Mediflo Duo is for single patient use. To prevent the potential transmission of infection, do not share Mediflo Duo with others.
- Do not use Mediflo Duo simultaneously with other physical activities.
- Do not use if there are any unintended sharp edges on the product.
- Hazards can result from the unauthorised modification to device. Do not modify this equipment without authorization from the manufacturer.

Cleaning Instructions

- After each use, rinse the device with a mild soap solution at lukewarm temperature and let air-dry.
- Once the components are dry, reassemble them and store in resealable plastic bag. The use of cleaning solutions which are not recommended by Flexicare may result in damage to the device and harm to patient.

Storage

Store in a cool, dry place out of direct sunlight.

Disposal

When the device shows signs of damage/degradation, dispose of in accordance with local policies and / or as advised by your local healthcare provider.

Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State of the user and/or patient.

Inspiration training according to the SMI method

SMI (Sustained Maximal Inspiration - in other words, breathing in slowly, continuously and as deeply as possible) is a conscious inhalation breathing exercise which corresponds to sighing and yawning at a physiological level.

Clinical benefits

The SMI method (slow continuous and deep inhalation) encourages even gas distribution in patients whose breathing resistance varies in different parts of the lungs and provides for faster recovery of the pre-operative status of surgical patients suffering from major restrictions in ventilation following thoracotomies or epigastric laparotomies. The SMI method has proved its worth as a physiological breathing exercise to prevent atelectasis and pneumonia for patients both in the surgical and internal medicine sector.

How is Mediflo Duo prepared and adjusted for SMI exercises?

Mediflo Duo has been factory-set for inspiration exercises (SMI), i.e. it is ready for immediate usage.

- White plug is already inserted in vertical stub.
- Arrow points to "Inspiration" (blue scale).
- Corrugated tube is already fitted to horizontal stub.

To adjust the Mediflo Duo to the ability of the patient: turn the cover on the blue scale to select the value between 1 and 5 that gives the patient maximum inhalation for 3 - 5 seconds. During this process the indicator ball should float continuously without striking the cover.

Performance of SMI exercises

1. We greatly recommend familiarizing the patient with Mediflo Duo at the pre-operative stage. The patient's upper body should be as upright as possible (set up back support). One hand holds the Mediflo Duo and the other the mouthpiece. The patient breathes in through the mouthpiece slowly, continuously and as deeply as possible so that the ball floats for 3 - 5 seconds (= SMI manoeuvre). The patient should hold his breath briefly, remove

the mouthpiece from his mouth and breathe out without making any effort. The patient should then breathe in and out calmly several times bypassing the Mediflo Duo.

2. The patient then takes a long deep breath in through the mouthpiece and breathes out again by passing the device.

One training session of 10 SMI manoeuvres is recommended per hour. Many therapists ask the patient to perform 10 SMI manoeuvres one after the other. With the provision of detailed instruction, inspiration exercises frequently, by themselves and of their own accord.

At the post-operative stage the scale value is initially reduced by around half to prevent the patient from becoming discouraged or overtiring himself. The scale value is increased as training progresses.

Expiration training according to the PEP method

PEP (Positive Expiratory Pressure) is an active (but not forced) manoeuvre involving breathing out against a resistance.

Clinical benefits

Breathing out against a resistance increases intrapulmonary pressure and not only reduces tracheobronchial collapse with unstable bronchial systems but also improves the transport of secretions. This method is therefore often used for patients with COPD to mobilise and eliminate secretions and prevent atelectasis.

How is Mediflo Duo prepared and adjusted for PEP exercises?

Mediflo Duo can be set up as a PEP training device in four steps

- Remove corrugated tube from horizontal stub.
- Remove white plug from vertical stub and insert in horizontal stub.
- Turn metering cover so arrow points to "Expiration" (orange scale).
- Fit corrugated tube to vertical stub.

To adjust Mediflo Duo to the ability of the patient: turn the cover to select the value on the orange scale that allows the patient to carry out the exercises for 2 minutes without undue effort. The higher the value, the shorter expiration time will be for the same expiration volume. The floating indicator ball provides an expiratory pressure of approx. 5cmH2O.

Performance of PEP exercises

1. The patient's upper body should be as upright as possible (set up back support). One hand holds the Mediflo Duo and the other the mouthpiece.
2. After a deep breath the patient breathes out into the unit, thus making the indicator ball float (expiratory pressure approx. 5cmH2O).
3. He breathes in either via the unit or with an open mouth and/or through the nose, bypassing the unit.
4. If he also breathes in through the mouthpiece, a negative inspiratory pressure of approx. 5 cm H2O is produced. This supports the ventilation of peripheral lung areas. This exercise is repeated for up to 5 minutes, several times a day.

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**Flexicare Medical Limited**
Flexicare Deutschland GmbH
Flexicare Medical Switzerland GmbH
c/o BK-Services AG Suurstoffi 22, 6343