Y-Piece Breathing Trainer


CE UK
CA

EN	Y-Piece Breathing Trainer
DE	Y-Atemtrainer
FR	Entraîneur respiratoire en Y
IT	Allenatore della respirazione Y-Piece
NL	Y-stuk ademtrainer
KO	Y-피스 호흡 훈련기

PEP

SMI

White Cap (Level 1)

Yellow Cap (Level 2)

Green Cap (Level 3)

Red Cap (Level 4)

105mm

EN – Y-PIECE BREATHING TRAINER

INTENDED PURPOSE

Breathing Trainers are manually operated incentive spirometers intended for inhalation training according to the Sustained Maximal Inspiration (SMI) method and expiration training according to the Positive Expiratory Pressure (PEP) method.

They are intended for use in adult and paediatric patients by medical professionals, carers and patients in hospital and out-of-hospital settings. Breathing Trainers are non-sterile single patient use devices, to be used for up to 30 days in a hospital environment and up to 6 months in a home environment.

INDICATIONS

Indications for the use of breathing trainers include:

- Acute and chronic respiratory failure caused by surgery.
- Neurological or musculoskeletal dysfunction, e.g., MND.
- Increasing lung volume by increasing FRC and Vt.
- Reducing hyperinflation / air trapping in e.g.: emphysema, bronchitis, asthma.
- Improving airway clearance in; cystic fibrosis (over 4-year-old); chronic bronchitis; bronchiectasis; bronchiolitis obliterans.
- Maximising the delivery of bronchodilators in patients receiving bronchial hygiene therapy.

CONTRAINDICATIONS

Contraindications for the use of breathing trainers include:

- Untreated pneumothorax
- Intracranial pressure > 20mm Hg
- Active haemoptysis
- Recent trauma or surgery to the skull, face, mouth, or oesophagus
- Patient with acute asthma attack or acute worsening of Chronic Obstructive Pulmonary Disease (COPD), unable to tolerate increased work of breathing
- Acute sinusitis or epistaxis
- Tympanic membrane rupture or other known or suspected inner ear pathology
- Nausea
- Following coronary artery bypass graft or valve replacement surgery
- When intended user:
 - is unable to generate adequate inspiration with a vital capacity < 10 mL/kg or an inspiratory capacity < 33% of predicted normal range.

EN – Y-PIECE BREATHING TRAINER

- unable to deep breathe effectively due to pain, diaphragmatic dysfunction, or opiate analgesia
- cannot be instructed or supervised to assure appropriate use of the device, or where co-operation and user understanding of the device is absent.
- is confused, delirious, heavily sedated or comatose.

CLINICAL BENEFITS

The SMI method (slow continuous and deep inhalation) encourages even gas distribution in patients whose breathing resistance varies in different parts of the lungs and provides for faster recovery of the pre-operative status of surgical patients suffering from major restrictions in ventilation following thoracotomies or epigastric laparotomies.

The SMI method has proved its worth as a physiological breathing exercise to prevent atelectasis and pneumonia for patients both in the surgical and internal medicine sector.

The PEP method (breathing out against a resistance) increases intrapulmonary pressure and not only reduces tracheobronchial collapse with unstable bronchial systems but also improves the transport of secretions. This method aids patients with COPD and other diseases to mobilise and eliminate secretions and prevent atelectasis.

WARNINGS

- Do not use if there is damage, discoloration or anything that may reduce the functionality of the device.
- Do not use if package has been previously opened or damaged.
- Do not use after expiry date.
- Do not stop taking any medication unless advised by a healthcare professional.
- Do not use Y-Piece Breathing Trainer during other physical activities.
- Do not modify this equipment without authorization from the manufacturer. Hazards can result from the unauthorised modification to device.
- Do not share Y-Piece Breathing Trainer with others to prevent the potential transmission of infection. Device may be used multiple times but may only be used by a single patient.
- Conduct a visual inspection on all of the components before use and do not use if the product is damaged.

EN – Y-PIECE BREATHING TRAINER**INSTRUCTIONS FOR USE****BEFORE USE**

- Consult a healthcare professional before use if the patient is on medication or suffers from a medical condition.
- The patient should be familiarised with Y-Piece Breathing Trainer before use.
- This Y-Piece Breathing Trainer has been developed for both inhalation and exhalation therapy.
 - Perform a visual inspection to ensure that all valves are properly seated. Failure to confirm correct valve positioning may lead to device not functioning as intended.
- The inspiratory limb is labelled with an arrow to indicate flow direction.
- The performance enhancement is achieved using an artificial brake.
- Breath training can be divided into four different performance ranges, for which there are four stenosis caps with increasing levels of artificial brake: 1st level = white cap, 2nd level = yellow cap, 3rd level = green cap, 4th level = red cap.

DURING USE**SMI METHOD (SUSTAINED MAXIMAL INSPIRATION)**

1. Select a suitable stenosis cap as recommended by a healthcare professional and place it on the inspiratory limb socket.
2. The patient's upper body should be as upright as possible (set up back support).
3. The patient breathes in through the mouthpiece slowly, continuously, and as deeply as possible (= SMI manoeuvre).
4. The patient should hold their breath briefly, remove the mouthpiece from their mouth and breathe out without making any effort.
5. The patient should then calmly breathe in and out several times bypassing the Y-Piece Breathing Trainer.
6. Repeat steps 3-5. The number of SMI manoeuvres should be dictated by a healthcare professional.

PEP METHOD (POSITIVE EXPIRATORY PRESSURE)

1. Select a suitable stenosis cap as recommended by a healthcare professional and place it on the expiratory limb socket.
2. The patient's upper body should be as upright as possible (set up back support).
3. The patient breathes out through the mouthpiece slowly, continuously, and as deeply as possible.
4. The patient should hold their breath briefly, remove the mouthpiece from their mouth and breathe in without making any effort.

EN – Y-PIECE BREATHING TRAINER

5. The patient should then calmly breathe in and out several times by passing the Y-Piece Breathing Trainer.
6. Repeat steps 3-5. The number of PEP manoeuvres should be dictated by a healthcare professional.

REUSE (CLEANING INSTRUCTIONS)

After each use, rinse the device with a mild soap solution at lukewarm temperature and let air-dry or disinfect using 70% iso-propyl alcohol solution and let air-dry.

Once the components are dry, store in resealable bag. The use of cleaning solutions which are not recommended by Flexicare may result in damage to the device and harm to patient.

AFTER USE (DISPOSAL)

When the device shows signs of damage/degradation, dispose of in accordance with local policies, and/or as advised by your local healthcare provider.

WHEN TO CONSULT A HEALTHCARE PROFESSIONAL

Any serious incident that occurs in relation to the device should be reported to the manufacturer, and the competent authority of the Member State of the user and/or patient.

STORAGE AND TRANSPORT

Transport and store in a cool, dry place out of direct sunlight.

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**Flexicare Medical Limited**
**Flexicare EU Ops Ltd.**
**Flexicare Medical Switzerland GmbH**
c/o BK-Services AG Suurstoffi 22